

# Treadmills

Clinical & Sports Medicine

Full range perfectly adapted for use in Clinical and Sports medicine

**Athlete Model 870 C**

**Sport Model 870 S**

**Clinical Model 870 A**



## Key Point

The running belt is triple thickness and is designed to absorb shock making for full comfort to the exercise.

## Features

- Parallel bars for the subject safety in event of fatigue or fall.
- Re-enforces the rigidity of the unit, increased mechanical strength.

## Options

- PC Console with 11" Touch Screen
- Belt for HR Monitoring with Bluetooth Technology
- PC ECG Medcard
- Saturometer SAO2
- Lanyard Auto-Stop
- Short Hand Rail (Only 870A model)
- Structure support for safety harness with auto-stop
- Paediatric Hand Rail

## Structure supports for safety harness with auto-stop



## PC Console



## Paediatric Hand Rail



## Short Hand Rail



## Technical Specifications

Models	Clinical 870 A	Sport 870 S	Athlete 870 C
Width	205 - 80 cm	203 - 89 cm	252 - 96 cm
Height	112 cm	115 cm	115 cm
Weight	163 Kg	220 Kg	245 Kg
Running Belt	50 - 142 cm	58 - 170 cm	65 - 190 cm
Height from floor	22 cm	16 cm	16 cm
Permissible Load	200 Kg	200 Kg	180 Kg
Speed Range	0 - 20 Km/h	0 - 25 Km/h	0 - 32 Km/h
Speed Increments	0.1 Km/h	0.1 Km/h	1 Km/h
Incline Range	0 - 25 %	0 - 22 % (25% opt)	0 - 22 % (25% opt)
Incline Increments	0.5 %	0.5 %	0.5 %
Speed Motor	Asynchronous 3 phases 2.2 HP, AC	Asynchronous 3 phases 2.2 HP, AC	Asynchronous 3 phases 2.2 HP, AC
Elevation Motor	90 Watt DC	90 Watt DC	90 Watt DC
Voltage Supply	230 V 50 HZ (opt 115 V AC, 60 HZ)	230 V 50 HZ (opt 115 V AC, 60 HZ)	230 V 50 HZ (opt 115 V AC, 60 HZ)
Max Continuous Power	1.5 KVA	1.5 KVA	1.5 KVA
Peak Power	5.0 KVA	5.0 KVA	5.0 KVA
Computer Interface	RS 232 or USB	RS 232 or USB	RS 232 or USB
Electrical Safety	Compliant w. EC 60601-1	Compliant w. EC 60601-1	Compliant w. EC 60601-1

This product is manufactured in accordance with EN-ISO 13485 (2012)

### Love Medical Ltd

Unit 1, Willan Enterprise Centre,  
Fourth Avenue, Village, Trafford Park,  
Manchester, M17 1DB, United Kingdom

[www.lovemedical.com](http://www.lovemedical.com)  
[info@lovemedical.com](mailto:info@lovemedical.com)  
[sales@lovemedical.com](mailto:sales@lovemedical.com)

Tel: 0161 976 2744  
Fax: 0161 976 2764